## **TERMINAL UNIT - MEDIUNIT (CZ)**





Recessed version



Exposed version



Installation plug

Medical terminal units provide quick and easy connection of hospital ward gas equipment to the hospital gas source. The type of medical gasoutlets are decided by national standards in each country and sometimes from local requests in each hospital. GCE complies with ISO 7396 and national installation standards with secure products where every product is fully tested in production. Our Medical gas outlets are in accordance with ISO EN 9170-1, ISO EN 9170-2 international standards.

**GCC** druva MED

- Wall housing is compatible with all GCE MediUnit standards like DIN, BSI, SS, CZ
- All functional components are from brass
- Simple installation
- Fast connection and disconnection
- Designed for medical environment, Small size and Easy to clean
- Complies with colour coding and description by standard
- After 10 years it is possible to upgrade the units with a special upgrade pack
- Recessed, Exposed and Bed head installation versions

Item No.	Denomination	Туре	Marking
0732030	O <sub>2</sub> – RECESSED	Pipe ø 10 mm	02
0732031	AIR – RECESSED	Pipe ø 10 mm	AIR
0732032	VAC – RECESSED	Pipe ø 10 mm	VAC
0732033	N <sub>2</sub> O – RECESSED	Pipe ø 10 mm	N <sub>2</sub> O
0732034	CO <sub>2</sub> – RECESSED	Pipe ø 10 mm	CO <sub>2</sub>
0732035	O <sub>2</sub> – EXPOSED	Pipe ø 10 mm	O <sub>2</sub>
0732036	AIR – EXPOSED	Pipe ø 10 mm	AIR
0732037	VAC – EXPOSED	Pipe ø 10 mm	VAC
0732038	N <sub>2</sub> O – EXPOSED	Pipe ø 10 mm	N <sub>2</sub> O
0732039	CO <sub>2</sub> – EXPOSED	Pipe ø 10 mm	CO <sub>2</sub>


## **INSTALLATION TOOLS**

Item No.	Denomination	
MP_00345	QC installation keys	
MP_00324	Button remover	
MP_01157ST	Pendants/bedhead unit – installation tool	
0732040	Installation plug (10 pcs)	

## TECHNICAL DATA

Gases:	O <sub>2</sub> ; N <sub>2</sub> O; Air, CO <sub>2</sub> ; AGSS; VAC
Dimensions:	Height: 73 mm, Width: 73 mm, Depth: 63 mm
Working pressure:	4–5 bar (breathing gases)
	7–10 bar (instrumental gases)
	(-0,4)–(-0,9) bar (vacuum)
Maximum pressure:	20 bar
Regulatory status:	Complies with Medical Devices Directive 93/42/EEC
	Complies with EN ISO 7396-1 (Central Gas Supply Systems)
	Complies with EN ISO 9170-1 (Terminal units)
	Complies with EN ISO 9170-2 (Terminal units for AGSS)
	Complies with CSN 85 2762 (Czech gas specific connections)
	present HTM 02-01

## **BASIC DIMENSIONS**



Note! Measurements in mm.